

Gusta™ & Gusta™ NC Product Specification:

Usable length	1410mm	Balloon coating	Hydrophilic(With Selective Coating Technology)
Proximal shaft	1.8F	Balloon material	Pebax/Nylon
Proximal shaft material	Stainless Steel	Fold	Tri-fold or Qua-fold
Distal shaft	2.7F	Guidewire	Maximum wire diameter 0.014"
Distal shaft material	Pebax/Nylon	Crossing profile	0.025"
Balloon markers	2 (1 or 2 markers for 1.5mm and 2.0mm)	Tip lead in profile	0.017" (0.0432mm)

Gusta™ Gusta™ NC

balloon dilatation catheter

Gusta™ Ordering Information

Balloon diameter(mm)	Balloon length(mm)				
	10	15	20	25	30
1.50	DBC1510	DBC1515	DBC1520	DBC1525	DBC1530
2.00	DBC2010	DBC2015	DBC2020	DBC2025	DBC2030
2.25	DBC2210	DBC2215	DBC2220	DBC2225	DBC2230
2.50	DBC2510	DBC2515	DBC2520	DBC2525	DBC2530
2.75	DBC2710	DBC2715	DBC2720	DBC2725	DBC2730
3.00	DBC3010	DBC3015	DBC3020	DBC3025	DBC3030
3.25	DBC3210	DBC3215	DBC3220	DBC3225	DBC3230
3.50	DBC3510	DBC3515	DBC3520	DBC3525	DBC3530

Gusta™ NC Ordering Information

Balloon diameter(mm)	Balloon length(mm)					
	8	10	12	15	18	20
2.25	DNC2208	DNC2210	DNC2212	DNC2215	DNC2218	DNC2220
2.50	DNC2508	DNC2510	DNC2512	DNC2515	DNC2518	DNC2520
2.75	DNC2708	DNC2710	DNC2712	DNC2715	DNC2718	DNC2720
3.00	DNC3008	DNC3010	DNC3012	DNC3015	DNC3018	DNC3020
3.25	DNC3208	DNC3210	DNC3212	DNC3215	DNC3218	DNC3220
3.50	DNC3508	DNC3510	DNC3512	DNC3515	DNC3518	DNC3520
3.75	DNC3708	DNC3710	DNC3712	DNC3715	DNC3718	DNC3720
4.00	DNC4008	DNC4010	DNC4012	DNC4015	DNC4018	DNC4020
4.25	DNC4208	DNC4210	DNC4212	DNC4215	DNC4218	DNC4220
4.50	DNC4508	DNC4510	DNC4512	DNC4515	DNC4518	DNC4520



Email: Marketing@demaxmedical.com


Beijing Office
Tel: +86 (10) 67868057
Fax: +86 (10) 67860415
Address: 1st Floor, Building No.1, No.8 Hongdabei road, BDA, Beijing China

www.demaxmedical.com

Sydney Office
Tel: +61 (02) 9939 7093
Fax: +61 (02) 9939 0397
Address: Suite 241, 117 Old Pittwater Rd Brookvale, NSW, 2100




With Hearts™




Gusta™

- Easy lesion access
- Sepcial designed tip forming




- Precise balloon expand compliance guarantee the exact dilatation
- Extremely low profile for CTO
- Enhanced welding technology for crossability



®

Gusta™ NC

- Suitable for 6F kissing balloon technique



* Test data on file at Demax Medical.

- Retentive pleating and folding
- Effective lubricious coating for maneuverability through repeated crossings
- Perfect for multi-dilatation



External Coating

- Selective Coating Technology