



Inflation Device

Selok[®] technology inside

Unique ergonomical design

Consistant pressure output

Precise manometer



CE₀₁₂₃ ARTG SFDA

Ref No.	Description
DID30	Max output pressure 30atm, 20ml volume
DID30S(ON)	Inflation Device, Max output pressure 30atm, 20ml volume Including High Pressure 3-way stopcock, "ON" Handle
DID30S(OFF)	Inflation Device, Max output pressure 30atm, 20ml volume, Including High Pressure 3-way stopcock, "OFF" Handle

International Marketing

Tel: +61 2 800 617 46

Fax: +61 2 808 864 77

Email: marketing@demaxmedical.com